

PERCUSSOR™

Hand-Held Percussion System



ERCHONIA
World Leaders in Low Level Laser Technology



TECHNICAL INFORMATION:

Variable Speed 6 to 60 percussions per second.

Patented Galvanic Motion

Heat Dispensing Electronics

Padded Velcro Grip

CNC Machined Medical Grade Delrin Plastic

Made in an ISO Certified Facility in the USA

100-120 Vac, 50/60Hz only

Summary:

Trigger-Point Therapy

Scar Tissue Mobilization

Comes with 2 attachments, Flat Pad and Cone Head & Carrying Case

Roller attachments & Trigger Point Series sold separately

Two Year Warranty

Compliant with:

SOR 98/282 Part 1, Sections 10-21

ISO 14971

EMC 2004/108/EC

LVD 2006/95/EC

IEC 60601-1-2 EMC

The ERCHONIA® PERCUSSOR is a patented industry leader in hand held percussion systems. With variable speeds up to 60 percussions a second and a heavy duty power package, it's treatment versatility can be tailored to the need of your vast patient base. The Percussor is manufactured in the USA with the highest quality materials in an ISO certified facility, thus the 2-year warranty we endorse on this proven therapy instrument. This professional and unique adjunct tool comes complete and ready to use right out of the box with no assembly required; simply select the attachment and plug into a power outlet. *Revitalize your practice with the most superior percussion therapy system instantly.*

The Percussor in Action!



Trigger Point Pad



Spinal Roller



Trigger Point Series



VP Roller



For more info, Call toll free: 1-888-242-0571
or visit www.erchonia.com